EN  
BEN  
DE  
BEN  
NL  
BEN  
FR  
BEN  
ES  
SON  
PER  
IT •  
BEN  
NO  
BEN  
SV  
BEN  
DA  
BEN  
FI •  
BEN  
PL  
BEN  
TR  
BEN  
EL •  
BEN  
CS  
BEN  
SK  
BEN  
HU  
BEN  
HR  
BEN  
ET  
BEN

**Manufacturer:** Nutricia Medical Devices BV, Taurusavenue 167, 2132 LS Hoofddorp, The Netherlands, [www.nutriciaflocares.com](http://www.nutriciaflocares.com)  
**Production location:** Nutricia Pharmaceutical Wuxi Co Ltd., P.R. China

**Man**  
**Pro**

CE 0344 **STERILE EO MD** 200345\_C June 2022

03

outer cover back page (1)

outer

NUTRICIA

flocare®

## FLOCARE® BENGMARK® PEG/J

### EN • INSTRUCTIONS FOR USE

BENGMARK® PEG/J

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BENGMARK® PEG/J

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BENGMARK® PEG/J

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BENGMARK® PEG/J

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


outer cover back page (2)

outer


outer cover right page to be folded inside (3)

inner cover page, to be folded inside: (4)

inner c

GBR

IRL

AUS

NZL

DEU

AUT

CHE

NLD

BEL

FRA

ESP

ITA

NOR

SWE

DNK

FIN

POL

TUR

GRC

CZE

SVK

HUN

HRV

EST

WWW

inner cover page middle (5) opposite of outer cover front page (2)

Inner

|            |   |
|------------|---|
| <b>GBR</b> | Nutricia Ltd, White Horse Business Park, TROWBRIDGE, Wiltshire, BA14 0XQ, England. Tel: 01225 711 688, <a href="http://www.nutricia.co.uk">www.nutricia.co.uk</a> Patient/Carer Information Number (UK only) Tel: 0800 093 3672 |
| <b>IRL</b> | Nutricia Ireland Ltd., Block 1, Deansgrange Business Park, Deansgrange, Co. DUBLIN, Ireland. <a href="http://www.nutricia.ie">www.nutricia.ie</a>   |
| <b>AUS</b> | Nutricia Australia Pty Ltd, Lvl 4, Bld D, 12-24 Talavera Rd, Macquarie Park NSW 2113, Australia, 1800 060 051, <a href="http://www.nutriciamedical.com.au">www.nutriciamedical.com.au</a>                                       |
| <b>NZL</b> | Danone Nutricia NZ Ltd, 56-58 Aintree Avenue, Airport Oaks, Auckland 2022, New Zealand, 0800 636 228, <a href="http://www.nutriciamedical.co.nz">www.nutriciamedical.co.nz</a>  |
| <b>DEU</b> | Nutricia Milupa GmbH, D - 60329 Frankfurt/Main, Deutschland, <a href="http://www.nutricia.de">www.nutricia.de</a><br>Kontakt Nutricia Milupa D-A-CH: 00800-70050000;  |
| <b>AUT</b> | Nutricia Milupa GmbH, A - 5412 Puch bei Hallein, Österreich, <a href="http://www.nutricia.at">www.nutricia.at</a><br>Kontakt Nutricia Milupa D-A-CH: 00800-70050000;  |
| <b>CHE</b> |  Nutricia Milupa SA, 8005 Zürich, Schweiz, <a href="http://www.nutricia.ch">www.nutricia.ch</a> Kontakt Nutricia Milupa D-A-CH: 00800-70050000 |
| <b>NLD</b> | Nutricia Nederland B.V., Antwoordnummer 10230, 2700 VB Zoetermeer-NL, Nutricia Medische Voedingsservice: 0800-022 33 22 (gratis).   |
| <b>BEL</b> | N.V. Nutricia België/Nutricia Belgique S.A., Werkhuizenkaai/Quai des Usines 160 - 1000 Brussel/Bruxelles, België/Belgique; 0800 99 486 (gratis/gratuit) - <a href="http://www.nutricia.be">www.nutricia.be</a>                  |
| <b>FRA</b> | NUTRICIA Nutrition Clinique, 17/19 rue des deux gares, CS 50149, 92565 Rueil-Malmaison Cedex, France. <a href="http://www.nutricia.fr">www.nutricia.fr</a>  |
| <b>ESP</b> | Distribuido en España por: Danone Nutricia S.R.L., C/ Torrelaguna 77, 6ª plta. CP 28043, Madrid. Tel. Atención al cliente: (+34) 900 211 088. <a href="http://www.nutricia.es">www.nutricia.es</a>                              |
| <b>ITA</b> | Distribuito da: Danone Nutricia S.p.A. Società Benefit, via C. Farini 41, 20159 Milano. Per informazioni: 800-822096 <a href="http://www.nutricia.it">www.nutricia.it</a>   |
| <b>NOR</b> | Nutricia Norge AS, Drammensveien 123, 0277 OSLO   |
| <b>SWE</b> | Nutricia Nordica AB, Gårdsvägen 14, 169 70 Solna  |
| <b>DNK</b> | Nutricia A/S, Rørmosevej 2A, 3450 Allerød   |
| <b>FIN</b> | Nutricia Medical Oy, Linnankatu 26 A, 20100 TURKU   |
| <b>POL</b> | Nutricia Polska Sp. z o.o., ul. Bobrowiecka 8, 00-728 Warszawa  |
| <b>TUR</b> | Numil Gıda Ürünleri San. ve Tic. A.. Yap Kredi Plaza B Blok No:1B Kat:8-9 Cömert Sok. Levent Mah. Beikta-STANBUL  |
| <b>GRG</b> | Numil Hellas A.E. 17 χλμ. Αθηνών Λαμίας & Καλαμπάς 2, 14564 Κηφισιά Γραμμή φροντίδας:800 11 68600, 210 6248 547   |
| <b>CZE</b> | Danone a.s., V parku 2294/2, 148 00 Praha 4 - Chodov, infolinka: 800 110 001, <a href="http://www.nutricia.cz">www.nutricia.cz</a>  |
| <b>SVK</b> | Danone s.r.o., Mlynské nivy 5, 821 09 Bratislava - mestská časť Ružinov. Infolinka: 0800 444 006  |
| <b>HUN</b> | Forgalmazza: Danone Kft., 1134 Budapest, Váci út 35. <a href="http://www.nutricia.hu">www.nutricia.hu</a>   |
| <b>HRV</b> | Distributer za Hrvatsku: G-M Pharma Zagreb d.o.o., Velika cesta 74, 10020 ZAGREB, Republika Hrvatska, Tel: 01/6234-053, <a href="http://www.g-m-pharma.hr">www.g-m-pharma.hr</a>  |
| <b>EST</b> | NutriMedical OÜ, TALLINN 11313, Eesti   |

[www.nutriciaflocare.com](http://www.nutriciaflocare.com)

Inner cover page (6) (opposite of outer cover page (1))

# EN • FLOCARE® BENCHMARK® PEG/J

## PRODUCT DESCRIPTION

The Flocare® Benchmark® PEG/J is a jejunal feeding tube which should be used with a Flocare® PEG CH 18. Nutrition/medication can be administered using an ENFit enteral giving set or enteral syringe.

## CONTENTS

|                               |  |
|-------------------------------|--|
| Tube material:                | blue fully opaque polyurethane tubing (A)  |
| External diameter:            | CH 9   |
| Length:                       | 105cm  |
| Tip:                          | antenna-tip with 2 lateral outlets (B).  |
| Benchmark®-spiral (patented): | 2,5 loops with a diameter of approx. 3 cm and a length of approx. 23 cm. These loops are straightened with a guide wire during insertion. After removal of the guide wire, the spiral allows optimal stable positioning in the small intestine. A black square shape marking printed on the tube about 25 cm from the tip indicates to which point the spiral will roll up when the guide wire is removed.   |
| Guide wire:                   | a prelubricated metallic guide wire, partially inserted into the tube (C).   |
| ENFit Connector:              | consisting of lilac (D) and turquoise (E) parts, to enable connection of the jejunal tube with the Y-connector.  |
| Y-Connector (red):            | to allow simultaneous aspiration of gastric contents and feeding in the small intestine. The gastric leg (F) contains a funnel connector allowing connection of a step connector with a drainage bag. The jejunal leg (G) is pre-assembled on the jejunal tube. The Y-connector contains an inner ring (H) to attach the connector to the PEG, as well as a separate red ring (I) to fixate the connection. The turquoise ring (J) on top of the Y-connector enables connection with the jejunal tube. |
| Tuohy borst adapter:          | to ensure optimal connection after tube placement; consisting of white, transparent (K) and turquoise (L) components which can be disconnected separately.   |
| Distance marks:               | every 10 centimeters.  |

## INDICATIONS

The Flocare® Benchmark® PEG/J CH 9 in combination with the Flocare® PEG CH 18 is suitable for those patients requiring enteral tube feeding directly in the jejunum or duodenum via a gastrostomy with a jejunal feeding tube, potentially together with simultaneous gastric drainage.

Indications for postpyloric tube feeding:

- Surgical patients
- Burn patients
- Patients with increased risk of aspiration or reflux and vomiting (early postoperative phase)
- Patients with limited gastric motility
- Patients with acute pancreatitis

This group of patients typically have a functioning intestinal tract with impaired stomach motility and/or increased risk of aspiration, e.g. early postoperative feeding. Early postoperative intestinal tube feeding reduces the risk of atrophy of the small intestine and reduces the risk of bacterial translocation.

It is recommended to deliver tube feed using an enteral feeding pump, to control the flow rate.

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For postpyloric feeding:

- Intestinal absorption failure
- Paralytic ileus
- Acute abdomen

More detailed information can be found in the Instructions For Use for the Flocare® PEG.

Notice: any serious incident that has occurred in relation to this device is to be reported to Nutricia and the local competent authority.

## PLACEMENT INSTRUCTIONS

PEG and PEG/J placement is only to be executed by trained and qualified healthcare professionals.

Preparation of the Flocare • PEG CH 18

After initial placement of the Flocare® PEG CH 18, with the endoscope still in stomach:

- Cut-off PEG tube at 35-cm marking.
- Pull external retention disc over the tube, but do not fix the tube in the disc yet.
- Place quick release clamp on the tube.
- Lubricate the inner-lumen of the PEG tube with approx. 10 ml water, using a syringe.

### A. Endoscopic placement:

1. Insert the guide wire completely into the jejunal tube. Withdraw the guide wire about 10-20 cm, and then push the guide wire until it reaches the end of the tube. Attach the handle of the guide wire (C) firmly to the ENFit Connector (D).
2. Lubricate the stretched jejunal tube with water to facilitate intubation.
3. Push the antenna tip (B) of the jejunal tube through the PEG CH 18 until the tip is visible in the stomach. If necessary add some extra water through the PEG CH 18 tube.
4. Grasp the antenna-tip with the endoscope forceps (fig. 1) and with help of the endoscope guide the tube through the pylorus. The tube tip should be positioned as far down the small intestine as possible, preferably beyond the Treitz ligament (fig. 2).
5. Keep the tube tip in place with the forceps, whilst withdrawing the endoscope (withdrawal of the endoscope approximately 30 cm is possible) (fig. 3).
6. Whilst the tube is still held by the forceps, withdraw the guide wire about 25 cm from the tube. Now the spiral will begin to form. With the forceps still in place, move the blue jejunal tube forward further into the small intestine (fig. 4).
7. Finally remove the endoscope, the forceps and the guide wire cautiously (fig. 5) being careful not to move the tube backwards at the same time.
8. Detach the ENFit connector from the Y-connector by unscrewing the turquoise ring (J) on top of the Y-connector (fig. 6).
9. Unscrew the separate red ring (I) from the Y-connector and move it over the PEG CH18 (fig. 7).
10. Connect the Y-connector to the PEG CH18 by fitting the red inner shaft (H) tightly into the PEG tube. Screw the separate red ring (I) into the Y-connector to fixate the Y-connector on the PEG tube (fig. 8)
11. Cut the ENFit connector from the jejunal tube (5 cm above the Y-connector) (fig. 9). If necessary, check the position of the tip of the tube endoscopically, to rule out any dislocation at this time.
12. Unscrew the turquoise ring (L) of the Tuohy borst adapter from the transparent part (K). Slide the turquoise ring over the blue jejunal tube. Fix the jejunal tube over the metallic pin inside the transparent part of the Tuohy borst adapter. Make sure the tube is completely moved over the metallic pin. Connect the turquoise and transparent parts of the Tuohy borst adapter by screwing in the turquoise part. Afterwards, connect the complete assembled connector to the Y-connector, by rotating the turquoise ring of the Y-connector (J). (fig. 10 and fig. 11).

13. Fix the PEG CH 18 tubing in the external fixation disc. Ensure the correct positioning of the external fixation disc is maintained by placing the blue safety clamp immediately behind the disc.
14. Check that the jejunal tube works by flushing with approx. 20 ml water.
15. Check if the jejunal tube is in the correct position with an X-ray. Note the length of the remaining part of the intestinal tube.

#### B. Fluoroscopic placement:

Follow point 1-3 of the instructions mentioned under A.

4. Check that the tube tip is in the stomach with fluoroscopy and direct the tip towards the pylorus. Withdraw the guide wire approximately 2-3 cm which will cause the tube to bend. This can facilitate finding the pylorus and help inserting the tube-tip through it. The safety tip design of the Bengmark® PEG/J tube enables manipulation with the guide wire in situ without the risk of lateral guide wire penetration, which could cause tissue damage. After pylorus penetration the guide wire should be brought back in the original position.
5. Insert the Bengmark® PEG/J tube with the guide wire as far as possible in the small intestine, preferably beyond the Treitz ligament.
6. Once the tube is in place, withdraw the guide wire.
7. Follow point 8-15 of the instructions mentioned under A.

### MAINTENANCE AND AFTER CARE

Verify status of the tube by checking the position and state of the external fixation disc and the tube itself, each time you change the tube feed, if in doubt whether the tube is in the right position, and at least three times a day. The Flocare® Bengmark® PEG/ J is used solely in combination with the Flocare® PEG CH 18, which in turn implies similar maintenance practices. Ensure you follow the maintenance instructions for the PEG carefully, also when used in combination with the Bengmark® PEG/ J. Be aware that maintenance with the PEG/J - PEG combination differs in one area (see below).

Within the PEG instructions for use the following is mentioned: Release the tube out of the external fixation disc, turn the tube daily 180° on its axis, and move tube up and down in the fully healed stoma (min 15 mm). When the PEG is combined with a Bengmark® PEG/J, the tube should NOT be turned on its axis, but only moved up and down in the fully healed stoma (min 15 mm).

Turning the tube is NOT carried out so as to avoid twisting of the Bengmark® PEG/J CH 9 tube, leading to possible dislocation of the PEG/J Bengmark® coil. Always return the tube to its original position with the PEG's blue safety clamp behind the external retention disk. (see figure 12).

Appropriate ancillaries, including enteral extension sets, giving sets and syringes, should be used using ENFit (enteral) connectors of a recognized standard (ISO80369-3) in order to avoid misconnection errors. Measures shall be taken to ensure correct drug dosing when administered through enteral tubes, for example using (low - dose tip) ENFit enteral syringes.

### WARNINGS

Disinfectants such as Povidone - Iodine (PVP -1; e.g.: Iso - Batadine®, Braunol®) and Octenidindihydrochlorid - Phenoxvethanol (e.g.: Octenisept®) should not be used because repeated exposure to this disinfectant might negatively affect the physical/mechanical properties of the tube. It is recommended to use disinfectants based on Polyhexanid (e.g. Prontosan®)

- The Flocare PEG/J Tuohy borst adaptor has a metal connection pin for fixation in the tube - do not expose the PEG/J to areas of strong magnetic fields, e.g. MRI equipment.
- There are no incompatibilities known in the choice of substances or (enteral) medicinal products administered via the PEG/J. Precaution is to be in place in delivering any kind of (enteral) medication via the PEG/J, as unknown effects - for other substances or higher concentrations - can occur. After any medication delivery, the PEG/J is to be flushed with water
- Material defects and tube occlusion can be largely ruled out if the Flocare® PEG/J is adequately flushed

before and after every administration of nutrition and/or medication

- Use liquid medication where possible
- Avoid using acidic substances such as cranberry juice and cola beverages to flush feeding tubes as the acidic quality when combined with formula proteins may actually contribute to tube clogging.

---

#### **GENERAL FLUSHING AND SYRINGE USE POLICY**

---

- Routinely flush the tube with water before and after the administration of feed or medication and at least every 8 hours to prevent blockage of the tube.
- As practicable, use an enteral syringe (minimum size of 20mL) for flushing the tube and always flush in a controlled manner.
- For medication administration directly into the tube, always administer using an appropriate size enteral syringe for the volume of medication but ensure delivery is in a controlled manner. In case more than one medication is administered at one time; make sure the tube is flushed thoroughly between types.
- Never use excessive force. If there is any resistance, stop the procedure and consult your clinician.

---

#### **DURATION OF USE**

---

The life span of the product varies according several factors including use of medication, gastric pH, patient motility and tube care.

The Tuohy borst adapter (K+L) is available as spare part and could extend the lifespan of the Flocare PEG/J or can replace a worn out or leaking connector within the normal foreseeable use period.

---

#### **SINGLE PATIENT USE**

---

Do not reuse this product - Do not reuse, reprocess or re-sterilize this device. Reuse, reprocessing or re-sterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness or death. Reuse, reprocessing or re-sterilization may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient.

---

#### **REMOVAL**

---

It is recommended to place a new Flocare® Benchmark® PEG/ J tube every 6 weeks.

Replacement or removal of the product may differ from this recommendation, based on medical scrutiny, viability of the product and/ or potential patient inconveniences and/or should be performed under medical instruction.

After use, dispose of product and packaging in accordance with hospital, administrative and/or local policy.

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