NUTRICIA  
**flocare**<sup>®</sup>

**FLOCARE<sup>®</sup> PEG SET**

**Manufacturer:** Nutricia Medical Devices BV, Taurusavenue 167, 2132 LS Hoofddorp, The Netherlands, <http://www.nutriciaflocare.com>  
**Production location:** Nutricia Pharmaceutical Wuxi Co Ltd., P.R. China

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
**ET • KASUTUSJUHEND**

PEG SET







**M** Manufacturer: Nutricia Medical Devices BV, Taurusavenue 167, 2132 LS Hoofddorp, The Netherlands, <http://www.nutriciaflocare.com>  
Production location: Nutricia Pharmaceutical Wuxi Co Ltd., P.R. China

CE  
0344



STERILE EO MD

200343\_D  
July 2022

outer cover right page to be folded inside (3)

<b>GBR</b>	Nutricia Ltd, White Horse Business Park, TROWBRIDGE, Wiltshire, BA14 0XQ, England. Tel: 01225 711 688, <a href="http://www.nutricia.co.uk">www.nutricia.co.uk</a> Patient/Carer Information Number (UK only) Tel: 0800 093 3672
<b>IRL</b>	Nutricia Ireland Ltd., Block 1, Deansgrange Business Park, Deansgrange, Co. DUBLIN, Ireland. <a href="http://www.nutricia.ie">www.nutricia.ie</a>
<b>AUS</b>	Nutricia Australia Pty Ltd, Lvl 4, Bld D, 12-24 Talavera Rd, Macquarie Park NSW 2113, Australia, 1800 060 051, <a href="http://www.nutriciamedical.com.au">www.nutriciamedical.com.au</a>
<b>NZL</b>	Danone Nutricia NZ Ltd, 56-58 Aintree Avenue, Airport Oaks, Auckland 2022, New Zealand, 0800 636 228, <a href="http://www.nutriciamedical.co.nz">www.nutriciamedical.co.nz</a>
<b>DEU</b>	Nutricia Milupa GmbH, D - 60329 Frankfurt/Main, Deutschland, <a href="http://www.nutricia.de">www.nutricia.de</a> Kontakt Nutricia Milupa D-A-CH: 00800-70050000;
<b>AUT</b>	Nutricia Milupa GmbH, A - 5412 Puch bei Hallein, Österreich, <a href="http://www.nutricia.at">www.nutricia.at</a> Kontakt Nutricia Milupa D-A-CH: 00800-70050000;
<b>CHE</b>	 Nutricia Milupa SA, 8005 Zürich, Schweiz, <a href="http://www.nutricia.ch">www.nutricia.ch</a> Kontakt Nutricia Milupa D-A-CH: 00800-70050000
<b>NLD</b>	Nutricia Nederland B.V., Antwoordnummer 10230, 2700 VB Zoetermeer-NL, Nutricia Medische Voedingsservice: 0800-022 33 22 (gratis).
<b>BEL</b>	N.V. Nutricia België/Nutricia Belgique S.A., Werkhuizenkaai/Quai des Usines 160 - 1000 Brussel/Bruxelles, België/Belgique; 0800 99 486 (gratis/gratuit) - <a href="http://www.nutricia.be">www.nutricia.be</a>
<b>FRA</b>	NUTRICIA Nutrition Clinique, 17/19 rue des deux gares, CS 50149, 92565 Rueil-Malmaison Cedex, France. <a href="http://www.nutricia.fr">www.nutricia.fr</a>
<b>ESP</b>	Distribuido en España por: Danone Nutricia S.R.L., C/ Torrelaguna 77, 6ª plta. CP 28043, Madrid. Tel. Atención al cliente: (+34) 900 211 088. <a href="http://www.nutricia.es">www.nutricia.es</a>
<b>ITA</b>	Distribuito da: Danone Nutricia S.p.A. Società Benefit, via C. Farini 41, 20159 Milano. Per informazioni: 800-822096 <a href="http://www.nutricia.it">www.nutricia.it</a>
<b>PRT</b>	Distribuído em Portugal por: Danone Nutricia Unipessoal, Lda.   Apartado 8106, EC Cabo Ruivo, 1802-001 Lisboa, Telef: 214 259 600   Linha Verde: 800 206 799   <a href="http://www.nutricia.pt">www.nutricia.pt</a>
<b>NOR</b>	Nutricia Norge AS, Drammensveien 123, 0277 OSLO
<b>SWE</b>	Nutricia Nordica AB, Gårdsvägen 14, 169 70 Solna
<b>DNK</b>	Nutricia A/S, Rørmosevej 2A, 3450 Allerød
<b>FIN</b>	Nutricia Medical Oy, Linnankatu 26 A, 20100 TURKU
<b>POL</b>	Nutricia Polska Sp. z o.o., ul. Bobrowiecka 8, 00-728 Warszawa
<b>TUR</b>	Numul Gıda Ürünleri San. ve Tic. A.. Yap Kredi Plaza B Blok No:1B Kat:8-9 Cömert Sok. Levent Mah. Beikta-STANBUL
<b>GRG</b>	Numul Hellas A.E. 17 χλμ. Αθηνών Λαμίας & Καλαμάτας 2, 14564 Κηφισιά Γραμμή φροντίδας:800 11 68600, 210 6248 547
<b>CZE</b>	Danone a.s., V parku 2294/2, 148 00 Praha 4 - Chodov, infolinka: 800 110 001, <a href="http://www.nutricia.cz">www.nutricia.cz</a>
<b>SVK</b>	Danone s.r.o., Mlynské nivy 5, 821 09 Bratislava - mestská časť Ružinov. Infolinka: 0800 444 006
<b>HUN</b>	Forgalmazza: Danone Kft., 1134 Budapest, Váci út 35. <a href="http://www.nutricia.hu">www.nutricia.hu</a>
<b>LVA</b>	NUTRICIA SIA, Gustava Zemgala gatve 74, Rga, LV-1039, Latvija
<b>EST</b>	NutriMedical OÜ, TALLINN 11313, Eesti

[www.nutriciaflocare.com](http://www.nutriciaflocare.com)

number

v

e

cia

s

,  
799

hert

μή

1,

inner cover page right (6) - opposite of back page 1.

# EN • FLOCARE® PEG SET

## PRODUCT DESCRIPTION

The Flocare PEG set is a percutaneous endoscopic gastrostomy set intended for patients requiring long term intragastric administration of tube feeding. The Flocare PEG is intended to be placed using the 'pull' technique, under endoscopic control. The PEG Summary of Safety and Clinical Performance (SSCP) can be found in EUDAMED at <https://ec.europa.eu/tools/eudamed>, or provided upon request

## CONTENTS

- A clear polyurethane tube (A), 40 cm in length, with
  - 3 radio - opaque lines
  - A three - winged silicone internal retention disc
  - A Hydromer® coated tip with loop, for easy passage through the abdominal wall
  - Cm markings up to 12 cm from the internal disc
- An external retention disc, made of silicone, for fixation and prevention of kinks in the tube, for maximum patient comfort (B)
- Blue safety clamp for required tube fixation (C)
- A scalpel (D)
- An introducer needle with a trocar (E)
- A wire with loop, for a convenient and secure connection with the tube loop (F)
- An ENFit feeding connector (G)
- A quick release clamp (H).
- A PEG Implant Card with instructions

The Flocare PEG tube is available in CH 10, CH 14 and CH 18, all are 40 cm in length. (B), (C), (G) and (H) are also available as spare parts

## INDICATIONS

- Long - term intragastric tube feeding (with gastric motility).
- Stomach decompression

## CONTRAINDICATIONS

- Mechanical obstruction of the GI tract (unless the tube placement is indicated for stomach decompression)
- Severe coagulation disorders
- Severe ascites &/or gastritis
- Ileus
- Anorexia nervosa or severe psychoses or advanced dementia
- Acute pancreatitis & peritonitis
- When it is impossible to bring the anterior part of the stomach against the abdominal wall (subtotal gastrectomy, ascites, hepatomegaly, etc.)
- Interposed organs, e.g. liver, colon
- Marked peritoneal carcinomatosis
- Extensive tumor infiltration in the area of the puncture site
- Oesophageal stenoses

Notice: any serious incident that has occurred in relation to this device is to be reported to Nutricia and the local competent authority

## PEG PLACEMENT

Preparation prior to endoscopic placement via Pull technique of Flocare PEG tube:

PEG tube placement is only to be executed by trained and qualified healthcare professionals.

The PEG tube insertion is performed with the patient in a supine position. The patient is prepared as usual for endoscopy. Thorough oral and pharyngeal hygiene is very important. Administration of a prophylactic antibiotic should be considered. If possible, explain the procedure to the patient.

In addition to the Flocare PEG set, the following equipment is required:

- Flexible endoscope with biopsy forceps
- Sterile gauze
- Skin disinfectant
- Syringe with local anaesthetic and long IM needle

PEG tube Placement:

1. Normally PEG tube placement is executed under mild sedation. If this is not the case, a topical anaesthetic of the vocal cords has to be considered. Pass the endoscope into the stomach and inflate the stomach. With the room lights dimmed, locate the appropriate puncture site in the left upper quadrant of the stomach, by placing a finger at the point of maximum transillumination by the endoscope light (fig.1). Palpate the abdominal wall to ensure the colon is not obstructing access to the stomach.
2. Clean the entire abdominal wall with a disinfectant and put sterile drapes around the puncture site. With the long IM needle, inject a local anaesthetic into the skin and in all layers of the abdominal wall. Insert the IM needle into the cavity of the stomach. Keep the IM needle in place by grasping it with the biopsy forceps. This avoids the stomach wall from dislocating versus the abdominal wall. (fig.2)
3. With the scalpel, make an incision of  $\pm 4 - 5$ mm in length very close to the IM needle. Push the puncture needle with the trocar through the incision into the stomach cavity (fig.3). Open the biopsy forceps, remove the IM needle and grasp the trocar with the forceps.
4. Remove the puncture needle and leave the trocar in place. Insert the pull wire through the trocar and grasp the wire with the biopsy forceps. (fig.4)
5. Dip the dilation tip of the tube in sterile water for approx. 10 seconds to activate the Hydromer® tipcoating. Pass the loop of the wire through the dilation loop, over the internal retention disc and pull snugly. (fig.5)
6. Pull the wire exiting the abdomen. The tube is then pulled down through the mouth into the stomach cavity and pulled out through the abdominal wall, using the dilation tip to dilate the exit site. The Hydromer® coating on the dilation tip will facilitate this. Remove the trocar when the conical tip of the tube is locked in the intragastric end of the trocar. Pull the internal retention disc of the tube snugly against the stomach wall (fig.6).
7. Cut the tube below the conical part, slide the silicone fixation disc over the tube and secure it by placing through the disc. Keep the external fixation disc tension - free against the skin for ca. 24h. Ensure the correct positioning of the external fixation disc is maintained by placing the blue safety clamp immediately behind the disc. Place the quick release clamp over the tube. Submerge the shaft of the feeding connector into water and push the shaft completely into the tube. The tube position should be confirmed endoscopically or by X - ray.
8. Record the brand name, the Charrière, the batch/LOT information and the introduced length of the PEG in to the patients' medical records, complete and hand over to the patient the enclosed PEG implant card

## MAINTENANCE AND AFTER CARE

---

### AFTERCARE POST ENDOSCOPIC PLACEMENT OF FLOCARE PEG SET:

- Tube feeding may commence 1 hour after successful PEG placement, however preferably 6 - 8 hours after PEG tube placement in order to prevent peritonitis
- Immediately after placement of the PEG, the external retention disc and safety clamp should be subjected to very low traction, without tension
- Allow a distance of approximately 5 mm (fig. 8) between the external retention disc and the skin, in order to place a sterile Y-compress under the external retention disc
- Inspect the stoma-site daily for any redness and swelling
- Disinfect the skin daily
- Cleansing and sterile renewal of dressings initially on a daily basis (usually up to day 7 post procedure and later 2-3 days). Keeping the stoma site dry is important to healing.

### MAINTENANCE OF FLOCARE PEG SET AFTER A HEALED STOMA SITE:

- As soon as the stoma is completely formed, it is sufficient to wash, rinse and dry the skin thoroughly. Inspect the condition of the PEG daily.
- Release the tube out of the external retention disc, turn the tube daily 180° on its axis, and move the tube up and down in the stoma (min 1.5 cm). After mobilization, return the tube to its initial position with some free distance (5 mm) between the skin and the external retention disc (fig.9).
- Always keep the blue safety clamp in place, also during wound care. The clamp indicates the correct position of the external retention disc and guarantees that the tube is fixed back in its original position after wound care.
- Check the position of the Flocare PEG by confirming the position of the external retention disc or, if possible, by measuring pH of aspirated gastric content. A pH value below 5.5 confirms correct position of the tube in the stomach. Do this each time you change the feed container, if there is any doubt concerning the position of the tube and at least 3 times a day.
- Appropriate ancillaries, including enteral extension sets, giving sets and syringes, should be used using ENFit (enteral) connectors of a recognized standard (ISO80369-3) in order to avoid misconnection errors.
- Measures shall be taken to ensure correct dosing when administered through enteral tubes, for example using (low - dose tip) ENFit enteral syringes

## GENERAL FLUSHING AND SYRINGE USE POLICY

---

- Routinely flush the tube with water before and after the administration of feed or medication and at least every 8 hours to prevent blockage of the tube.
- As practicable, use an enteral syringe (minimum size of 20mL) for flushing the tube and always flush in a controlled manner.
- For medication administration directly into the tube, always administer using an appropriate size enteral syringe for the volume of medication but ensure delivery is in a controlled manner. In case more than one medication is administered at one time; make sure the tube is flushed thoroughly between types.
- Never use excessive force. If there is any resistance, stop the procedure and consult your clinician

### WARNINGS:

---

- Any significant change to the PEG (discoloration, hairline cracks, leakage, etc.) should be reported to the attending physician.
- Always keep the blue safety clamp in place.

- Disinfectants such as Povidone - Iodine (PVP - I; e.g.: Iso - Betadine®, Braunol®) and Octenidindihydrochlorid - Phenoxyethanol (e.g.: Octenisept®) should not be used because repeated exposure to this disinfectant might negatively affect the physical/mechanical properties of the tube. It is recommended to use disinfectants based on Polyhexanid (e.g. Prontosan, Prontosan C.), based on Octenidine (e.g. Octenilin) or ActiMaris disinfectants wound irrigation solutions.
- The Flocare PEG connector CH10 has a metal connection pin for fixation in the PEG tube - do not expose the PEG CH10 with the CH10 ENFit connector attached to areas of strong magnetic fields, e.g. MRI equipment.
- There are no incompatibilities known in the choice of substances administered via the PEG. Precaution is to be in place in delivering via the PEG, as unknown effects can occur. After any administration delivery, the PEG is to be flushed with water
- Material defects and tube occlusion can be largely ruled out if the Flocare® PEG Set is adequately flushed before and after every administration of nutrition and/or medication
- Avoid using acidic substances such as cranberry juice and cola beverages to flush feeding tubes as the acidic quality when combined with formula proteins may actually contribute to tube clogging.
- The PEG has not been tested for use with medication. If alternative routes of drug delivery are not an option, the pharmacist or other knowledgeable healthcare professional needs to give authorization for medicines to be given in conjunction with enteral feeding and provide guidance on the preparation & administration to ensure safe use

#### DURATION OF USE

---

The PEG is designed to have a life span of several years, but the life span of the tube varies according several factors (see daily care and maintenance paragraph for more information). These factors include medication, gastric pH, patient motility and tube care. The durability of a PEG tube system is primarily linked to its careful handling. The following spare parts that could extend the lifespan of the Flocare PEG are available: External retention disk, quick release clamp, blue safety clamp and the ENFit feeding connector. The state and position of the tube should be checked endoscopically after 8 - 10 months

#### REMOVAL

---

- The tube should not be removed within 10 days of placement, or until the stoma tract is formed.
- Removal of this tube depends upon the viability of the product and/or should be performed under medical instruction.
- To remove the tube, cut it at skin level. The internal retention disc and the remaining tube must be retrieved by endoscope. After use, dispose of product and packaging in accordance with hospital, administrative and/or local policy.
- The Flocare PEG can be replaced by a replacement gastrostomy (Gastrostomy tube or Button).
- Oral feeding after PEG removal may only start after medical approval.

#### WARNINGS:

---

With patients who can no longer be subjected to endoscopy, or endoscopic removal method is not necessary, the PEG tube can be cut off at the skin level and removed, by allowing the internal retention disc to pass naturally. This does result in an increased risk of complications and bowel obstruction. Careful follow - up and monitoring of the patient is essential. This method is contraindicated in children and patients who have previously undergone surgery

## POTENTIAL COMPLICATIONS

---

Possible complications include, but may not be limited to:

Frequent potential complications:

- Local wound infection

Acute and severe potential complications that occur infrequently:

- Perforation
- Serious abdominal haemorrhage
- Peritonitis
- Aspiration
- Puncture of the colon
- Laceration of the gastric or small bowel mucosa, and stridor

Possible long - term complications:

- Tube occlusion
- Tube porosity and fracture with subsequent leakage from the tube or the tube connection
- Development of cellulitis
- Eczema or hypergranulation tissue

The development of most of these potential long - term complications is dependent on the quality of aftercare given to the tubing system and can be effectively avoided if the proper measures are taken.

Rare complication: buried bumper syndrome

## SINGLE PATIENT USE

---

Do not reuse, reprocess or re-sterilize this device. Reuse, reprocessing or re-sterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness or death. Reuse, reprocessing or re-sterilization may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient. After use, dispose of product and packaging in accordance with hospital, administrative and/or local government policy

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